

# User Manual:

## Endoscope Geometric Distortion evaluation software

### Tool Reference

RST Reference Number: RST24MD05.01

Date of Publication: 05/04/2023

Recommended Citation: U.S. Food and Drug Administration. (2023). Endoscope Geometric Distortion- A Method for Quantitative Evaluation (RST24MD05.01). <https://cdrh-rst.fda.gov/endoscope-geometric-distortion-method-quantitative-evaluation>

### For more information

[Catalog of Regulatory Science Tools to Help Assess New Medical Devices](#)

## Disclaimer

### About the Catalog of Regulatory Science Tools


The enclosed tool is part of the Catalog of Regulatory Science Tools, which provides a peer-reviewed resource for stakeholders to use where standards and qualified Medical Device Development Tools (MDDTs) do not yet exist. These tools do not replace FDA-recognized standards or MDDTs. This catalog collates a variety of regulatory science tools that the FDA's Center for Devices and Radiological Health's (CDRH) Office of Science and Engineering Labs (OSEL) developed. These tools use the most innovative science to support medical device development and patient access to safe and effective medical devices. If you are considering using a tool from this catalog in your marketing submissions, note that these tools have not been qualified as [Medical Device Development Tools](#) and the FDA has not evaluated the suitability of these tools within any specific context of use. You may [request feedback or meetings for medical device submissions](#) as part of the Q-Submission Program.


For more information about the Catalog of Regulatory Science Tools, email [OSEL\\_CDRH@fda.hhs.gov](mailto:OSEL_CDRH@fda.hhs.gov).

This manual provides instructions for utilizing the software designed to assess endoscope geometric distortion. MatLab (<https://www.mathworks.com/products/matlab.html>) is required to run this software. The software includes the following parts:


- A main program named “distortion\_polyfit\_correct\_7p5.m”.
- A folder named “private”. Several subroutines are saved in this folder and the main program will call them during running. *If you run the software on a Mac computer, this folder might not be recognized. In that case, you can move the subroutines to the same folder as the main program.*
- A folder named “images”. Images in this folder will be inputs to the software. A sample image named “OLYM0013.jpg” is in the folder for testing the software. This software can analyze images in jpg, tif, and png formats. You can put the target images taken by your own endoscopes into this folder as inputs. Please refer to our paper for taking target images. [1]
- A folder named “outputs”. An output file will be generated in this folder after each run. The file contains the values of all the variables for the analyzed endoscope image. To understand the algorithm and the meaning of these variables, please refer to our paper [1].

#### You can follow the following steps to run the software:



- Follow our paper [1] to take target images and put the images in the “images” folder.
- Install Matlab on your computer.
- Download the zip file named “GeometricDistortion\_7p5\_share\_2.zip” from GitHub and unzip it.
- Double click the main program named “distortion\_polyfit\_correct\_7p5.m” to open it in MatLab.
- Click the “Run” button  to run the software.
- The “images” folder will open. Select an image you want to analyze and click the “Open” button.



- A figure window will open. Follow the instruction on the top of the window to select cross points on the horizontal and vertical lines that across the image center as follows.
  - Starting at the center, click the left mouse button once at each intersection until the leftmost intersection is selected. Then press the Enter key on the keyboard to end the selection.
  - Repeat similar procedure from the center to right, from the center to top and from the center to bottom.



- The software will generate the following outputs.
  - A “\*.mat” file in the “outputs” folder. Double click the output file to open it with MatLab to check the values of all the variables from the image analysis. The file name is a string started with the name of the analyzed image followed by date, month, year, and time when the file was generated.
  - Two images in two MatLab image windows as follows. The left image shows the relation between distorted radius ( $r_d$ ) and undistorted radius ( $r_u$ ) in four directions: center to left, center to right, center to top and center to bottom. The average of these four directions is also displayed. The right image shows the local magnification curves in the radial and tangential directions. The values for these curves (**r\_d** and **MLR\_N** for the radial direction; **r\_d** and **MLT\_N** for the tangential direction) can be found in the “\*.mat” output file in the “outputs” folder.



#### Reference:

- [1] Q. Wang, W.-C. Cheng, N. Suresh, and H. Hua, "Development of the local magnification method for quantitative evaluation of endoscope geometric distortion," Journal of biomedical optics 21, 056003 (2016). doi: <https://doi.org/10.1117/1.JBO.21.5.056003>

#### Disclaimer:

*The mention of commercial products, their sources, or their use in connection with material reported herein is not to be construed as either an actual or implied endorsement of such products by the Department of Health and Human Services."*

\*For further information, please email [quanzeng.wang@fda.hhs.gov](mailto:quanzeng.wang@fda.hhs.gov) .